

PERCI-SAR System VERSION 3.31 UPGRADE

This offer is only for PERCI-SAR system users with version 3.21 or earlier!

- **Version 3.31 is a release that has a new waveform graphical format! Pressure – Flow Graphs are now visually aligned!**
(In all previous versions the Perci-SAR Retrieve mode, even though cursor measurements were aligned in the time domain, the waveform graphs were not visually aligned)
- **Version 3.31 creates a new AREA waveform plot and adds this plot to the graph display options in the Retrieve Mode. This plot is a computed area waveform for all pressure – flow data points.**
- **Version 3.31 now tests oral – nasal pressure relationships in VP function samples! This will help the user to determine when error conditions might cause invalid area measurements. (For example: when there is incomplete lip closure or a blocked oral catheter.)**
- **Version 3.31 – All manual documentation is now available via the new Help System. All manuals are also installed and available in PDF format.**



Ordering Info:

PERCI-SAR SYSTEM VERSION 3.31 UPGRADE

\$199.00

To: MicroTronics Corp.

88 VilCom Circle /Suite 165

Chapel Hill, NC 27514

(919) 929-2657

Fax: (919) 960-6864